

PRODUCT DATA INFORMATION

AIR-SAN

DESCRIPTION

A synergistic blend of novel disinfectants, alcohols and fragrances. An air sanitizer, disinfectant and deodorizing spray, used in hospitals, sick rooms, mortuaries, hotel foyers, airport waiting rooms, vehicles and wherever air borne germs and surface disinfecting is needed.

CHEMICAL PROPERTIES

BIODEGRADABLE	
pH	N/A
SG	0.79
Odour	Alcohol Lavender
Auto-ignition Temperature	432°C

VITAL INFORMATION

AIR-SAN is effective against many different viruses, bacteria and fungus, namely the following air borne and surface viruses: Pseudomonas aerogenosa SATCC Pse 16, Staphylococcus aureus SATCC Sta 10, Escherichia coli SATCC est 25, and many more.

When tested in accordance with SABS 1615 – 1995, the above organisms achieved a 99.9% kill within 5 minutes.

METHOD OF USE

For air sanitizing and deodorizing, stand in the middle of a room and spray towards the ceiling corners of the room. Only a quick burst is needed.
For surface disinfecting, spray a liberal amount on the area to be disinfected and allow the product to remain on the surface for at least 5 minutes.

SAFETY PRECAUTIONS

Contents are flammable and are under pressure. Keep in a cool place away from heat. If product comes into contact with...

SKIN	Wash with soap and water for 15min
EYES	Remove contact lenses. Flush with water or saline for 15min
INGESTION	Do not induce vomiting. If conscious give 1-2 glasses of water. Get medical attention

KEEP OUT OF REACH OF CHILDREN

FOR MORE INFORMATION CONTACT 082 873 2688